

Bellman & Symfon Deaf and Hard of Hearing Alarm System

Installation Directions

- 1) Select a location in the home to install the smoke alarm – See suggestions at bottom
- 2) Remove the mounting plate from the smoke alarm by twisting in counter clockwise direction. Remove back padding if applicable.
- 3) Place the mounting plate on the wall or ceiling and mark the location of the two screw holes.
- 4) Use the supplied screws to attach mounting plate to wall or ceiling.
- 5) Connect the battery to the battery snaps on the back of the alarm and push down into slot.
 - a. Red and green lights on the front of the alarm should begin to blink every 45 seconds
- 6) Place the smoke alarm on top of the baseplate and turn clockwise to attach
- 7) Press and hold test button to test the alarm
- 8) Get out the flash receiver
- 9) Plug flash receiver into wall outlet in bedroom (near the bed)
- 10) Pull the plastic strip at the bottom of the flash alarm to turn on
- 11) Press and hold the test button on the flash alarm until green and yellow lights blink
- 12) Press and hold the test button on the smoke alarm until it sounds
- 13) The flash receiver will indicate that pairing was successful by switching between LED lights 4 times.
- 14) Get out the bed shaker
- 15) Plug bed shaker into the back of the flash receiver
- 16) Press and hold test button on the smoke alarm to confirm that pairing of all devices was successful. (Smoke alarm should sound, flash receiver should flash, and bed shaker should shake)
- 17) Place bed shaker underneath pillow or mattress
- 18) Instruct resident to test alarm system at least once a month

Note:

Flash receiver and bed shaker will turn off after 10 seconds when being tested

Smoke alarm location suggestions:

It is recommended to install the smoke alarm outside the main sleeping area, in a hallway or corridor. Install alarm at least 12 inches away from corners, light fixtures, and air ducts.

